

Floor-mounted Analog X-ray System

YJ-FA100U

Technical Parameters

Model	YJ-FA100U-32	YJ-FA100U-50	YJ-FA100U-65
Power Output	32kW	50kW	65kW
Current	10-400mA	10-630mA	10-800mA
Voltage	40kV-150kV	40kV-150kV	40kV-150kV
mAs	0.1-400mAs	0.1-630mAs	0.1-800mAs
Exposure time	1-6300ms	1-6300ms	1-6300ms
Focus Size	0.6/1.2	0.6/1.2	0.6/1.2



Floor-mounted Analog X-ray System adopts floor-mounted table and wall stands. Its smart design, flexible movement, easy operation and high working efficiency provide more imaging positions.

Gallery



C-spine



Chest



L-spine



Pelvis



Hands



Foot



Knee joint



L-spine



Chest



Tube Head

Specially designed tube head arm holder makes smooth vertical movement and compensates mechanical vibration.



Tube Stand

Saving the cost and installation time designed by freestanding (Floor mounted). Accommodates facilities where can't support a ceiling mounted system, Wide range of travels covering lateral and other difficult positioning. Electromagnetic locks secure at any position.



Wall Bucky Stand

Wide vertical-travel range with minimal focal spot to floor distance for comfortable standing position. Electromagnetic locks secure at any vertical position. Effortless unforced motions by counterbalanced. Cassette flexibility with durable bucky tray.

